



Pwortho[®]



Product
Catalogue

HANGZHOU PENGWU MEDICAL EQUIPMENT CO.,LTD

PENGWU ORTHODONTICS CO.,LTD (HONG KONG)

E-mail:alfred@pw-ortho.com T:+86 13916823991

[Http://www.pw-ortho.com](http://www.pw-ortho.com)



Company Introduction

PengWu Orthodontics is a manufacturer of high quality orthodontic products, manufacturing and selling self-ligating metal brackets, conventional metal brackets, buccal tubes, band, arch wires, elastomer and other accessories, all products achieved CE, FDA, ISO13485 certificate.

With over 10 years experience in orthodontic industry, PengWu Orthodontics has set up complete MIM CIM CNC production lines that make our products to be at the advanced level in China and all over the world. Our products have been used in many famous hospitals and exported to over 30 countries.

PengWu Orthodontics has established good cooperation relationships with well-known universities and hospitals for better design, products innovation, such as Beijing University, China Jiliang University, Shanghai Ninth People Hospital etc.

With high quality products and good services, PengWu would like to continue to build and keep long term business relationship and partnerships with customers around the world.

Orthodontic Brackets -----04

Buccal Tubes -----08

Molar Band -----12

Molar Band Prewelded with Tubes-13

Orthodontic Accessories-----15

Passive Self-Ligating Metal Brackets

- Clips works smoothly, easy open and close
- MIM technology, 17-4 stainless steel, precise and low friction force slot
- Low profile and low tie wings design with smooth surface, make patients more comfortable
- Laser welding, reliable and secure
- 80 gauge mesh compound contoured base provides a precise fit
- More efficient by reducing practice time
- Speed up treatment process
- Roth and MBT prescriptions are available
- improved oral health by elimination of elastic ties



Active Self-Ligating Metal Brackets

- Spring clip easy open and close
- Precise CNC technology, 17-4 stainless steel, precise and low friction force slot
- Laser welding, reliable and secure
- 80 gauge mesh compound contoured base provides a precise fit.
- Low profile and low tie wings design with smooth surface, make patients more comfortable
- Passive Phase, Expressive Phase, Active Phase greatly shorten the treatment time
- Roth and MBT prescriptions are available
- Improved oral health by elimination of elastic ties



Regular Metal Brackets

- MIM technology and CNC technology options are available
- Low profile design with rounded edge and smooth surface, enhance patient's comfort
- 80 gauge mesh compound contoured base provides a precise fit.
- Built in torque and angulations for accurate treating force transaction
- Vertical scribe line for easy positioning during practice
- Color ID for easy identification
- Mesh base laser marking option is available



Edgewise Brackets

Product Description	Hooks	Slot	Item Number (MIM)	Item Number (CNC)
Standard Edgewise Bracket	No Hook	0.018	BR131-20A	BR031-20A
		0.022	BR132-20A	BR032-20A
	3 with Hook	0.018	BR131-23A	BR031-23A
		0.022	BR132-23A	BR032-23A
	345 with Hook	0.018	BR131-25A	BR031-25A
		0.022	BR132-25A	BR032-25A
Mini Edgewise Bracket	No Hook	0.018	BR131-10A	BR031-10A
		0.022	BR132-10A	BR032-10A
	3 with Hook	0.018	BR131-13A	BR031-13A
		0.022	BR132-13A	BR032-13A
	345 with Hook	0.018	BR131-15A	BR031-15A
		0.022	BR132-15A	BR032-15A

Maxillary	Standard									
	UR5	UR4	UR3	UR2	UR1	UL1	UL2	UL3	UL4	UL5
Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Tip	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
In/Out	0.7	0.7	0.7	0.7	0.7	0.7	0.7	0.7	0.7	0.7
Width	3.3	3.3	3.6	3	3.3	3.3	3	3.6	3	3
Mandibular	Standard									
	LR5	LR4	LR3	LR2	LR1	LL1	LL2	LL3	LL4	LL5
Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Tip	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
In/Out	0.7	0.7	0.7	0.7	0.7	0.7	0.7	0.7	0.7	0.7
Width	3.3	3.3	3.6	2.5	2.5	2.5	2.5	3.6	3.3	3.3
Maxillary	Mini									
	UR5	UR4	UR3	UR2	UR1	UL1	UL2	UL3	UL4	UL5
Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Tip	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
In/Out	0.7	0.7	0.7	0.7	0.7	0.7	0.7	0.7	0.7	0.7
Width	3.2	3.2	3.3	2.7	3	3	2.7	3.3	3.2	3.2
Mandibular	Mini									
	LR5	LR4	LR3	LR2	LR1	LL1	LL2	LL3	LL4	LL5
Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Tip	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
In/Out	0.7	0.7	0.7	0.7	0.7	0.7	0.7	0.7	0.7	0.7
Width	3.2	3.2	3.3	2.3	2.3	2.3	2.3	3.3	3.2	3.2

Bracket Identification

upper right	upper left
lower right	lower left

No marking for Lower Anteriors



Roth Brackets

Product Description	Hooks	Slot	Item Number (MIM)	Item Number (CNC)
Standard Roth Bracket	No Hook	0.018	BR111-20A	BR011-20A
		0.022	BR112-20A	BR012-20A
	3 with Hook	0.018	BR111-23A	BR011-23A
		0.022	BR112-23A	BR012-23A
	345 with Hook	0.018	BR111-25A	BR011-25A
		0.022	BR112-25A	BR012-25A
Mini Roth Bracket	No Hook	0.018	BR111-10A	BR011-10A
		0.022	BR112-10A	BR012-10A
	3 with Hook	0.018	BR111-13A	BR011-13A
		0.022	BR112-13A	BR012-13A
	345 with Hook	0.018	BR111-15A	BR011-15A
		0.022	BR112-15A	BR012-15A
Passive Self-Ligating Roth Bracket	No Hook	0.018	BP111-10A	/
		0.022	BP112-10A	/
	3 with Hook	0.018	BP111-13A	/
		0.022	BP112-13A	/
	345 with Hook	0.018	BP111-15A	/
		0.022	BP112-15A	/
Active Self-Ligating Roth Bracket	No Hook	0.018	/	BA211-10A
		0.022	/	BA212-10A
	3 with Hook	0.018	/	BA211-13A
		0.022	/	BA212-13A
	345 with Hook	0.018	/	BA211-15A
		0.022	/	BA212-15A

Maxillary	Standard									
	UR5	UR4	UR3	UR2	UR1	UL1	UL2	UL3	UL4	UL5
Torque	-7°	-7°	-2°	8°	12°	12°	8°	-2°	-7°	-7°
Tip	0°	0°	11°	9°	5°	5°	9°	11°	0°	0°
In/Out	0.8	0.8	0.8	1.18	0.8	0.8	1.18	0.8	0.8	0.8
Width	3.2	3.2	3.3	2.9	3.6	3.6	2.9	3.3	3.2	3.2

Mandibular	Standard									
	LR5	LR4	LR3	LR2	LR1	LL1	LL2	LL3	LL4	LL5
Torque	-22°	-17°	-11°	-1°	-1°	-1°	-1°	-11°	-17°	-22°
Tip	0°	0°	5°	0°	0°	0°	0°	5°	0°	0°
In/Out	0.8	0.8	0.8	1.18	1.18	1.18	1.18	0.8	0.8	0.8
Width	3.2	3.2	3.2	2.5	2.5	2.5	2.5	3.2	3.2	3.2

Maxillary	Mini									
	UR5	UR4	UR3	UR2	UR1	UL1	UL2	UL3	UL4	UL5
Torque	-7°	-7°	-2°	8°	12°	12°	8°	-2°	-7°	-7°
Tip	0°	0°	11°	9°	5°	5°	9°	11°	0°	0°
In/Out	0.68	0.68	0.66	1.06	0.68	0.68	1.06	0.66	0.68	0.68
Width	2.8	2.8	3.2	2.8	3.2	3.2	2.8	3.2	2.8	2.8

Mandibular	Mini									
	LR5	LR4	LR3	LR2	LR1	LL1	LL2	LL3	LL4	LL5
Torque	-22°	-17°	-11°	-1°	-1°	-1°	-1°	-11°	-17°	-22°
Tip	0°	0°	5°	0°	0°	0°	0°	5°	0°	0°
In/Out	0.68	0.72	0.68	1.06	1.06	1.06	1.06	0.68	0.72	0.68
Width	2.8	2.8	2.8	2.5	2.5	2.5	2.5	2.8	2.8	2.8

Bracket Identification

upper right upper left
lower right lower left

No marking for Lower Anteriors

MBT Brackets

Product Description	Hooks	Slot	Item Number (MIM)	Item Number (CNC)
Standard MBT Bracket	No Hook	0.018	BR121-20A	BR021-20A
		0.022	BR122-20A	BR022-20A
	3 with Hook	0.018	BR121-23A	BR021-23A
		0.022	BR122-23A	BR022-23A
	345 with Hook	0.018	BR121-25A	BR021-25A
		0.022	BR122-25A	BR022-25A
Mini MBT Bracket	No Hook	0.018	BR121-10A	BR021-10A
		0.022	BR122-10A	BR022-10A
	3 with Hook	0.018	BR121-13A	BR021-13A
		0.022	BR122-13A	BR022-13A
	345 with Hook	0.018	BR121-15A	BR021-15A
		0.022	BR122-15A	BR022-15A
Passive Self-Ligating MBT Bracket	No Hook	0.018	BP121-10A	/
		0.022	BP122-10A	/
	3 with Hook	0.018	BP121-13A	/
		0.022	BP122-13A	/
	345 with Hook	0.018	BP121-15A	/
		0.022	BP122-15A	/
Active Self-Ligating MBT Bracket	No Hook	0.018	/	BA221-10A
		0.022	/	BA222-10A
	3 with Hook	0.018	/	BA221-13A
		0.022	/	BA222-13A
	345 with Hook	0.018	/	BA221-15A
		0.022	/	BA222-15A

Maxillary	Standard									
	UR5	UR4	UR3	UR2	UR1	UL1	UL2	UL3	UL4	UL5
Torque	-7°	-7°	0°	10°	17°	17°	10°	0°	-7°	-7°
Tip	0°	0°	8°	8°	4°	4°	8°	8°	0°	0°
In/Out	0.9	0.7	0.7	1.1	0.7	0.7	1.1	0.7	0.7	0.9
Width	3.2	3.2	3.3	2.9	3.6	3.6	2.9	3.3	3.2	3.2

Mandibular	Standard									
	LR5	LR4	LR3	LR2	LR1	LL1	LL2	LL3	LL4	LL5
Torque	-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
Tip	2°	2°	3°	0°	0°	0°	0°	3°	2°	2°
In/Out	0.7	0.7	0.7	1.1	1.1	1.1	1.1	0.7	0.7	0.7
Width	3.2	3.2	3.2	2.7	2.7	2.7	2.7	3.2	3.2	3.2

Maxillary	Mini									
	UR5	UR4	UR3	UR2	UR1	UL1	UL2	UL3	UL4	UL5
Torque	-7°	-7°	0°	10°	17°	17°	10°	0°	-7°	-7°
Tip	0°	0°	8°	8°	4°	4°	8°	8°	0°	0°
In/Out	0.7	0.7	0.7	1.08	0.7	0.7	1.08	0.7	0.7	0.7
Width	2.8	2.8	2.8	2.8	3.1	3.1	2.8	2.8	2.8	2.8

Mandibular	Mini									
	LR5	LR4	LR3	LR2	LR1	LL1	LL2	LL3	LL4	LL5
Torque	-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
Tip	2°	2°	3°	0°	0°	0°	0°	3°	2°	2°
In/Out	0.7	0.7	0.7	0.9	0.9	0.9	0.9	0.7	0.7	0.7
Width	2.8	2.8	2.8	2.2	2.2	2.2	2.2	2.8	2.8	2.8

Bracket Identification

upper right	upper left
lower right	lower left

No marking for Lower Anteriors

Buccal Tubes

- Low profile design and smooth surface, round corner enhances patient comfort
- Parelled smooth and round hook avoid gingival irritation and improve patient comfort
- Funnelled entracne helps to guide the wire into the buccal tube
- Round funnel entrance and square funnel entrance are available
- Fully contoured mesh base anatomically conforms to the teeth
- Torque built in design provide level slot line-up and help avoid occlusal interference
- 80 gauge mesh base ensures maximum bonding strength
- Sepcial notch on occlusal side make it easy to hold and place during practice
- Pronounced groove increases bonding accuracy and perfect fit with the tooth

1st Molar Buccal Tubes Single

50 pcs / pack (individual)

Prescription	Teeth	Torque	Offset	Width (mm)	Bondable (Round Funnel Entrance / Square Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6°	4	TU231-1126 / TU331-1126	TU232-1126 / TU332-1126	TU231-2126 / TU331-2126	TU232-2126 / TU332-2126
	UL6	0°	6°	4	TU231-1116 / TU331-1116	TU232-1116 / TU332-1116	TU231-2116 / TU331-2116	TU232-2116 / TU332-2116
	LL6	0°	6°	4	TU231-1136 / TU331-1136	TU232-1136 / TU332-1136	TU231-2136 / TU331-2136	TU232-2136 / TU332-2136
	LR6	0°	6°	4	TU231-1146 / TU331-1146	TU232-1146 / TU332-1146	TU231-2146 / TU331-2146	TU232-2146 / TU332-2146
Roth	UR6	-14°	10°	4	TU211-1126 / TU311-1126	TU212-1126 / TU312-1126	TU211-2126 / TU311-2126	TU212-2126 / TU312-2126
	UL6	-14°	10°	4	TU211-1116 / TU311-1116	TU212-1116 / TU312-1116	TU211-2116 / TU311-2116	TU212-2116 / TU312-2116
	LL6	-25°	4°	4	TU211-1136 / TU311-1136	TU212-1136 / TU312-1136	TU211-2136 / TU311-2136	TU212-2136 / TU312-2136
	LR6	-25°	4°	4	TU211-1146 / TU311-1146	TU212-1146 / TU312-1146	TU211-2146 / TU311-2146	TU212-2146 / TU312-2146
MBT	UR6	-14°	10°	4	TU221-1126 / TU321-1126	TU222-1126 / TU322-1126	TU221-2126 / TU321-2126	TU222-2126 / TU322-2126
	UL6	-14°	10°	4	TU221-1116 / TU321-1116	TU222-1116 / TU322-1116	TU221-2116 / TU321-2116	TU222-2116 / TU322-2116
	LL6	-20°	0°	4	TU221-1136 / TU321-1136	TU222-1136 / TU322-1136	TU221-2136 / TU321-2136	TU222-2136 / TU322-2136
	LR6	-20°	0°	4	TU221-1146 / TU321-1146	TU222-1146 / TU322-1146	TU221-2146 / TU321-2146	TU222-2146 / TU322-2146

50 pcs / pack (individual)

Prescription	Teeth	Torque	Offset	Width (mm)	Weldable (Regular: No Funnel Entracne / Round Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6°	4	TW131-1126 / TW231-1126	TW132-1126 / TW232-1126	TW131-2126 / TW231-2126	TW132-2126 / TW232-2126
	UL6	0°	6°	4	TW131-1116 / TW231-1116	TW132-1116 / TW232-1116	TW131-2116 / TW231-2116	TW132-2116 / TW232-2116
	LL6	0°	6°	4	TW131-1136 / TW231-1136	TW132-1136 / TW232-1136	TW131-2136 / TW231-2136	TW132-2136 / TW232-2136
	LR6	0°	6°	4	TW131-1146 / TW231-1146	TW132-1146 / TW232-1146	TW131-2146 / TW231-2146	TW132-2146 / TW232-2146
Roth	UR6	-14°	10°	4	TW111-1126 / TW211-1126	TW112-1126 / TW212-1126	TW111-2126 / TW211-2126	TW112-2126 / TW212-2126
	UL6	-14°	10°	4	TW111-1116 / TW211-1116	TW112-1116 / TW212-1116	TW111-2116 / TW211-2116	TW112-2116 / TW212-2116
	LL6	-25°	4°	4	TW111-1136 / TW211-1136	TW112-1136 / TW212-1136	TW111-2136 / TW211-2136	TW112-2136 / TW212-2136
	LR6	-25°	4°	4	TW111-1146 / TW211-1146	TW112-1146 / TW212-1146	TW111-2146 / TW211-2146	TW112-2146 / TW212-2146
MBT	UR6	-14°	10°	4	TW121-1126 / TW221-1126	TW122-1126 / TW222-1126	TW121-2126 / TW221-2126	TW122-2126 / TW222-2126
	UL6	-14°	10°	4	TW121-1116 / TW221-1116	TW122-1116 / TW222-1116	TW121-2116 / TW221-2116	TW122-2116 / TW222-2116
	LL6	-20°	0°	4	TW121-1136 / TW221-1136	TW122-1136 / TW222-1136	TW121-2136 / TW221-2136	TW122-2136 / TW222-2136
	LR6	-20°	0°	4	TW121-1146 / TW221-1146	TW122-1146 / TW222-1146	TW121-2146 / TW221-2146	TW122-2146 / TW222-2146

4 pcs / set

Prescription	Teeth	Torque	Offset	Width (mm)	Bondable (Round Funnel Entrance / Square Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6°	4	TU231-11B / TU331-11B	TU232-11B / TU332-11B	TU231-21B / TU331-21B	TU232-21B / TU332-21B
	UL6	0°	6°	4				
	LL6	0°	6°	4				
	LR6	0°	6°	4				
Roth	UR6	-14°	10°	4	TU211-11B / TU311-11B	TU212-11B / TU312-11B	TU211-21B / TU311-21B	TU212-21B / TU312-21B
	UL6	-14°	10°	4				
	LL6	-25°	4°	4				
	LR6	-25°	4°	4				
MBT	UR6	-14°	10°	4	TU221-11B / TU321-11B	TU222-11B / TU322-11B	TU221-21B / TU321-21B	TU222-21B / TU322-21B
	UL6	-14°	10°	4				
	LL6	-20°	0°	4				
	LR6	-20°	0°	4				

4 pcs / set

Prescription	Teeth	Torque	Offset	Width (mm)	Weldable (Regular: No Funnel Entracne / Round Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6°	4	TW131-11B / TW231-11B	TW132-11B / TW232-11B	TW131-21B / TW231-21B	TW132-21B / TW232-21B
	UL6	0°	6°	4				
	LL6	0°	6°	4				
	LR6	0°	6°	4				
Roth	UR6	-14°	10°	4	TW111-11B / TW211-11B	TW112-11B / TW212-11B	TW111-21B / TW211-21B	TW112-21B / TW212-21B
	UL6	-14°	10°	4				
	LL6	-25°	4°	4				
	LR6	-25°	4°	4				
MBT	UR6	-14°	10°	4	TW121-11B / TW221-11B	TW122-11B / TW222-11B	TW121-21B / TW221-21B	TW122-21B / TW222-21B
	UL6	-14°	10°	4				
	LL6	-20°	0°	4				
	LR6	-20°	0°	4				

1st Molar Buccal Tubes Double Combination Tubes

50 pcs / pack (individual)

Prescription	Teeth	Torque	Offset	Width (mm)	Bondable (Round Funnel Entrance / Square Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6°	4	TU231-1226 / TU331-1226	TU232-1226 / TU332-1226	TU231-2226 / TU331-2226	TU232-2226 / TU332-2226
	UL6	0°	6°	4	TU231-1216 / TU331-1216	TU232-1216 / TU332-1216	TU231-2216 / TU331-2216	TU232-2216 / TU332-2216
	LL6	0°	6°	4	TU231-1236 / TU331-1236	TU232-1236 / TU332-1236	TU231-2236 / TU331-2236	TU232-2236 / TU332-2236
	LR6	0°	6°	4	TU231-1246 / TU331-1246	TU232-1246 / TU332-1246	TU231-2246 / TU331-2246	TU232-2246 / TU332-2246
Roth	UR6	-14°	10°	4	TU211-1226 / TU311-1226	TU212-1226 / TU312-1226	TU211-2226 / TU311-2226	TU212-2226 / TU312-2226
	UL6	-14°	10°	4	TU211-1216 / TU311-1216	TU212-1216 / TU312-1216	TU211-2216 / TU311-2216	TU212-2216 / TU312-2216
MBT	LR6	-20°	0°	4	TU221-1226 / TU321-1226	TU222-1226 / TU322-1226	TU221-2226 / TU321-2226	TU222-2226 / TU322-2226
	LL6	-20°	0°	4	TU221-1216 / TU321-1216	TU222-1216 / TU322-1216	TU221-2216 / TU321-2216	TU222-2216 / TU322-2216

50 pcs / pack (individual)

Prescription	Teeth	Torque	Offset	Width (mm)	Weldable (Regular: No Funnel Entrance / Round Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6°	4	TW131-1226 / TW231-1226	TW132-1226 / TW232-1226	TW131-2226 / TW231-2226	TW132-2226 / TW232-2226
	UL6	0°	6°	4	TW131-1216 / TW231-1216	TW132-1216 / TW232-1216	TW131-2216 / TW231-2216	TW132-2216 / TW232-2216
	LL6	0°	6°	4	TW131-1236 / TW231-1236	TW132-1236 / TW232-1236	TW131-2236 / TW231-2236	TW132-2236 / TW232-2236
	LR6	0°	6°	4	TW131-1246 / TW231-1246	TW132-1246 / TW232-1246	TW131-2246 / TW231-2246	TW132-2246 / TW232-2246
Roth	UR6	-14°	10°	4	TW111-1226 / TW211-1226	TW112-1226 / TW212-1226	TW111-2226 / TW211-2226	TW112-2226 / TW212-2226
	UL6	-14°	10°	4	TW111-1216 / TW211-1216	TW112-1216 / TW212-1216	TW111-2216 / TW211-2216	TW112-2216 / TW212-2216
MBT	LR6	-20°	0°	4	TW121-1226 / TW221-1226	TW122-1226 / TW222-1226	TW121-2226 / TW221-2226	TW122-2226 / TW222-2226
	LL6	-20°	0°	4	TW121-1216 / TW221-1216	TW122-1216 / TW222-1216	TW121-2216 / TW221-2216	TW122-2216 / TW222-2216

4 pcs / set

Prescription	Teeth	Torque	Offset	Width (mm)	Bondable (Round Funnel Entrance / Square Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6°	4	TU231-12B / TU331-12B	TU232-12B / TU332-12B	TU231-22B / TU331-22B	TU232-22B / TU332-22B
	UL6	0°	6°	4				
	LL6	0°	6°	4				
	LR6	0°	6°	4				
Roth	UR6	-14°	10°	4	TU211-12B / TU311-12B	TU212-12B / TU312-12B	TU211-22B / TU311-22B	TU212-22B / TU312-22B
	UL6	-14°	10°	4				
	LL6	-25°	4°	4				
	LR6	-25°	4°	4				
MBT	UR6	-14°	10°	4	TU221-12B / TU321-12B	TU222-12B / TU322-12B	TU221-22B / TU321-22B	TU222-22B / TU322-22B
	UL6	-14°	10°	4				
	LL6	-20°	0°	4				
	LR6	-20°	0°	4				

4 pcs / set

Prescription	Teeth	Torque	Offset	Width (mm)	Weldable (Regular: No Funnel Entrance / Round Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6°	4	TW131-12B / TW231-12B	TW132-12B / TW232-12B	TW131-22B / TW231-22B	TW132-22B / TW232-22B
	UL6	0°	6°	4				
	LL6	0°	6°	4				
	LR6	0°	6°	4				
Roth	UR6	-14°	10°	4	TW111-12B / TW211-12B	TW112-12B / TW212-12B	TW111-22B / TW211-22B	TW112-22B / TW212-22B
	UL6	-14°	10°	4				
	LL6	-25°	4°	4				
	LR6	-25°	4°	4				
MBT	UR6	-14°	10°	4	TW121-12B / TW221-12B	TW122-12B / TW222-12B	TW121-22B / TW221-22B	TW122-22B / TW222-22B
	UL6	-14°	10°	4				
	LL6	-20°	0°	4				
	LR6	-20°	0°	4				



1st Molar Buccal Tubes Double Rectangular Tubes

50 pcs / pack (individual)

Prescription	Teeth	Torque	Offset	Width (mm)	Bondable (Round Funnel Entrance / Square Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6*	4	TU231-1226 / TU331-1226	TU232-1226 / TU332-1226	TU231-2226 / TU331-2226	TU232-2226 / TU332-2226
	UL6	0°	6*	4	TU231-1216 / TU331-1216	TU232-1216 / TU332-1216	TU231-2216 / TU331-2216	TU232-2216 / TU332-2216
	LL6	0°	6*	4	TU231-1236 / TU331-1236	TU232-1236 / TU332-1236	TU231-2236 / TU331-2236	TU232-2236 / TU332-2236
	LR6	0°	6*	4	TU231-1246 / TU331-1246	TU232-1246 / TU332-1246	TU231-2246 / TU331-2246	TU232-2246 / TU332-2246
Roth	UR6	-14°	10*	4	TU211-1226 / TU311-1226	TU212-1226 / TU312-1226	TU211-2226 / TU311-2226	TU212-2226 / TU312-2226
	UL6	-14°	10*	4	TU211-1216 / TU311-1216	TU212-1216 / TU312-1216	TU211-2216 / TU311-2216	TU212-2216 / TU312-2216
	LL6	-25°	4*	4	TU211-1236 / TU311-1236	TU212-1236 / TU312-1236	TU211-2236 / TU311-2236	TU212-2236 / TU312-2236
	LR6	-25°	4*	4	TU211-1246 / TU311-1246	TU212-1246 / TU312-1246	TU211-2246 / TU311-2246	TU212-2246 / TU312-2246
MBT	UR6	-14°	10*	4	TU221-1226 / TU321-1226	TU222-1226 / TU322-1226	TU221-2226 / TU321-2226	TU222-2226 / TU322-2226
	UL6	-14°	10*	4	TU221-1216 / TU321-1216	TU222-1216 / TU322-1216	TU221-2216 / TU321-2216	TU222-2216 / TU322-2216
	LL6	-20°	0*	4	TU221-1236 / TU321-1236	TU222-1236 / TU322-1236	TU221-2236 / TU321-2236	TU222-2236 / TU322-2236
	LR6	-20°	0*	4	TU221-1246 / TU321-1246	TU222-1246 / TU322-1246	TU221-2246 / TU321-2246	TU222-2246 / TU322-2246



50 pcs / pack (individual)

Prescription	Teeth	Torque	Offset	Width (mm)	Weldable (Regular: No Funnel Entrance / Round Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6*	4	TW131-1226 / TW231-1226	TW132-1226 / TW232-1226	TW131-2226 / TW231-2226	TW132-2226 / TW232-2226
	UL6	0°	6*	4	TW131-1216 / TW231-1216	TW132-1216 / TW232-1216	TW131-2216 / TW231-2216	TW132-2216 / TW232-2216
	LL6	0°	6*	4	TW131-1236 / TW231-1236	TW132-1236 / TW232-1236	TW131-2236 / TW231-2236	TW132-2236 / TW232-2236
	LR6	0°	6*	4	TW131-1246 / TW231-1246	TW132-1246 / TW232-1246	TW131-2246 / TW231-2246	TW132-2246 / TW232-2246
Roth	UR6	-14°	10*	4	TW111-1226 / TW211-1226	TW112-1226 / TW212-1226	TW111-2226 / TW211-2226	TW112-2226 / TW212-2226
	UL6	-14°	10*	4	TW111-1216 / TW211-1216	TW112-1216 / TW212-1216	TW111-2216 / TW211-2216	TW112-2216 / TW212-2216
	LL6	-25°	4*	4	TW111-1236 / TW211-1236	TW112-1236 / TW212-1236	TW111-2236 / TW211-2236	TW112-2236 / TW212-2236
	LR6	-25°	4*	4	TW111-1246 / TW211-1246	TW112-1246 / TW212-1246	TW111-2246 / TW211-2246	TW112-2246 / TW212-2246
MBT	UR6	-14°	10*	4	TW121-1226 / TW221-1226	TW122-1226 / TW222-1226	TW121-2226 / TW221-2226	TW122-2226 / TW222-2226
	UL6	-14°	10*	4	TW121-1216 / TW221-1216	TW122-1216 / TW222-1216	TW121-2216 / TW221-2216	TW122-2216 / TW222-2216
	LL6	-20°	0*	4	TW121-1236 / TW221-1236	TW122-1236 / TW222-1236	TW121-2236 / TW221-2236	TW122-2236 / TW222-2236
	LR6	-20°	0*	4	TW121-1246 / TW221-1246	TW122-1246 / TW222-1246	TW121-2246 / TW221-2246	TW122-2246 / TW222-2246



4 pcs / set

Prescription	Teeth	Torque	Offset	Width (mm)	Bondable (Round Funnel Entrance / Square Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6*	4	TU231-12B / TU331-12B	TU232-12B / TU332-12B	TU231-22B / TU331-22B	TU232-22B / TU332-22B
	UL6	0°	6*	4				
	LL6	0°	6*	4				
	LR6	0°	6*	4				
Roth	UR6	-14°	10*	4	TU211-12B / TU311-12B	TU212-12B / TU312-12B	TU211-22B / TU311-22B	TU212-22B / TU312-22B
	UL6	-14°	10*	4				
	LL6	-25°	4*	4				
	LR6	-25°	4*	4				
MBT	UR6	-14°	10*	4	TU221-12B / TU321-12B	TU222-12B / TU322-12B	TU221-22B / TU321-22B	TU222-22B / TU322-22B
	UL6	-14°	10*	4				
	LL6	-20°	0*	4				
	LR6	-20°	0*	4				



4 pcs / set

Prescription	Teeth	Torque	Offset	Width (mm)	Weldable (Regular: No Funnel Entrance / Round Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6*	4	TW131-12B / TW231-12B	TW132-12B / TW232-12B	TW131-22B / TW231-22B	TW132-22B / TW232-22B
	UL6	0°	6*	4				
	LL6	0°	6*	4				
	LR6	0°	6*	4				
Roth	UR6	-14°	10*	4	TW111-12B / TW211-12B	TW112-12B / TW212-12B	TW111-22B / TW211-22B	TW112-22B / TW212-22B
	UL6	-14°	10*	4				
	LL6	-25°	4*	4				
	LR6	-25°	4*	4				
MBT	UR6	-14°	10*	4	TW121-12B / TW221-12B	TW122-12B / TW222-12B	TW121-22B / TW221-22B	TW122-22B / TW222-22B
	UL6	-14°	10*	4				
	LL6	-20°	0*	4				
	LR6	-20°	0*	4				



50 pcs / pack (individual)

Prescription	Teeth	Torque	Offset	Width (mm)	Weldable (Regular: No Funnel Entrance / Round Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6*	4	TW131-1326 / TW231-1326	TW132-1326 / TW232-1326	TW131-2326 / TW231-2326	TW132-2326 / TW232-2326
	UL6	0°	6*	4	TW131-1316 / TW231-1316	TW132-1316 / TW232-1316	TW131-2316 / TW231-2316	TW132-2316 / TW232-2316
Roth	UR6	-14°	10*	4	TW111-1326 / TW211-1326	TW112-1326 / TW212-1326	TW111-2326 / TW211-2326	TW112-2326 / TW212-2326
	UL6	-14°	10*	4	TW111-1316 / TW211-1316	TW112-1316 / TW212-1316	TW111-2316 / TW211-2316	TW112-2316 / TW212-2316
MBT	LR6	-20°	0*	4	TW121-1326 / TW221-1326	TW122-1326 / TW222-1326	TW121-2326 / TW221-2326	TW122-2326 / TW222-2326
	LL6	-20°	0*	4	TW121-1316 / TW221-1316	TW122-1316 / TW222-1316	TW121-2316 / TW221-2316	TW122-2316 / TW222-2316



4 pcs / set

Prescription	Teeth	Torque	Offset	Width (mm)	Weldable (Regular: No Funnel Entrance / Round Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6*	4	TW131-13B / TW231-13B	TW132-13B / TW232-13B	TW131-23B / TW231-23B	TW132-23B / TW232-23B
	UL6	0°	6*	4				
Roth	UR6	-14°	10*	4	TW111-13B / TW211-13B	TW112-13B / TW212-13B	TW111-23B / TW211-23B	TW112-23B / TW212-23B
	UL6	-14°	10*	4				
MBT	LR6	-20°	0*	4	TW121-13B / TW221-13B	TW122-13B / TW222-13B	TW121-23B / TW221-23B	TW122-23B / TW222-23B
	LL6	-20°	0*	4				



2nd Molar Buccal Tubes Single

50 pcs / pack (individual) □

Prescription	Teeth	Torque	Offset	Width (mm)	Bondable (Round Funnel Entrance / Square Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egdewise	UR6	0°	6°	3.2	TU231-1127 / TU331-1127	TU232-1127 / TU332-1127	TU231-2127 / TU331-2127	TU232-2127 / TU332-2127
	UL6	0°	6°	3.2	TU231-1117 / TU331-1117	TU232-1117 / TU332-1117	TU231-2117 / TU331-2117	TU232-2117 / TU332-2117
	LL6	0°	6°	3.2	TU231-1137 / TU331-1137	TU232-1137 / TU332-1137	TU231-2137 / TU331-2137	TU232-2137 / TU332-2137
	LR6	0°	6°	3.2	TU231-1147 / TU331-1147	TU232-1147 / TU332-1147	TU231-2147 / TU331-2147	TU232-2147 / TU332-2147
Roth	UR6	-14°	10°	3.2	TU211-1127 / TU311-1127	TU212-1127 / TU312-1127	TU211-2127 / TU311-2127	TU212-2127 / TU312-2127
	UL6	-14°	10°	3.2	TU211-1117 / TU311-1117	TU212-1117 / TU312-1117	TU211-2117 / TU311-2117	TU212-2117 / TU312-2117
	LL6	-25°	4°	3.2	TU211-1137 / TU311-1137	TU212-1137 / TU312-1137	TU211-2137 / TU311-2137	TU212-2137 / TU312-2137
	LR6	-25°	4°	3.2	TU211-1147 / TU311-1147	TU212-1147 / TU312-1147	TU211-2147 / TU311-2147	TU212-2147 / TU312-2147
MBT	UR6	-14°	10°	3.2	TU221-1127 / TU321-1127	TU222-1127 / TU322-1127	TU221-2127 / TU321-2127	TU222-2127 / TU322-2127
	UL6	-14°	10°	3.2	TU221-1117 / TU321-1117	TU222-1117 / TU322-1117	TU221-2117 / TU321-2117	TU222-2117 / TU322-2117
	LL6	-10°	0°	3.2	TU221-1137 / TU321-1137	TU222-1137 / TU322-1137	TU221-2137 / TU321-2137	TU222-2137 / TU322-2137
	LR6	-10°	0°	3.2	TU221-1147 / TU321-1147	TU222-1147 / TU322-1147	TU221-2147 / TU321-2147	TU222-2147 / TU322-2147

50 pcs / pack (individual) □

Prescription	Teeth	Torque	Offset	Width (mm)	Weldable (Regular: No Funnel Entrance / Round Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egdewise	UR6	0°	6°	3.2	TW131-1127 / TW231-1127	TW132-1127 / TW232-1127	TW131-2127 / TW231-2127	TW132-2127 / TW232-2127
	UL6	0°	6°	3.2	TW131-1117 / TW231-1117	TW132-1117 / TW232-1117	TW131-2117 / TW231-2117	TW132-2117 / TW232-2117
	LL6	0°	6°	3.2	TW131-1137 / TW231-1137	TW132-1137 / TW232-1137	TW131-2137 / TW231-2137	TW132-2137 / TW232-2137
	LR6	0°	6°	3.2	TW131-1147 / TW231-1147	TW132-1147 / TW232-1147	TW131-2147 / TW231-2147	TW132-2147 / TW232-2147
Roth	UR6	-14°	10°	3.2	TW111-1127 / TW211-1127	TW112-1127 / TW212-1127	TW111-2127 / TW211-2127	TW112-2127 / TW212-2127
	UL6	-14°	10°	3.2	TW111-1117 / TW211-1117	TW112-1117 / TW212-1117	TW111-2117 / TW211-2117	TW112-2117 / TW212-2117
	LL6	-25°	4°	3.2	TW111-1137 / TW211-1137	TW112-1137 / TW212-1137	TW111-2137 / TW211-2137	TW112-2137 / TW212-2137
	LR6	-25°	4°	3.2	TW111-1147 / TW211-1147	TW112-1147 / TW212-1147	TW111-2147 / TW211-2147	TW112-2147 / TW212-2147
MBT	UR6	-14°	10°	3.2	TW121-1127 / TW221-1127	TW122-1127 / TW222-1127	TW121-2127 / TW221-2127	TW122-2127 / TW222-2127
	UL6	-14°	10°	3.2	TW121-1117 / TW221-1117	TW122-1117 / TW222-1117	TW121-2117 / TW221-2117	TW122-2117 / TW222-2117
	LL6	-10°	0°	3.2	TW121-1137 / TW221-1137	TW122-1137 / TW222-1137	TW121-2137 / TW221-2137	TW122-2137 / TW222-2137
	LR6	-10°	0°	3.2	TW121-1147 / TW221-1147	TW122-1147 / TW222-1147	TW121-2147 / TW221-2147	TW122-2147 / TW222-2147



Molar Band

Features

- Interior surface of the band is rough finished that ensures maximum bonding strength
- Anatomically contoured to tooth and edge hardness that provides excellent fit and stability
- Smooth surface and rolled edge enhances patient comfort
- Laser marks for easy permanent identification

Size	Perimeter (mm)		Item Number			
	Upper	Lower	UR	UL	LL	LR
28	28	28	BM6UR-280	BM6UL-280	BM6LL-280	BM6LR-280
28+	28.5	28.5	BM6UR-285	BM6UL-285	BM6LL-285	BM6LR-285
29	29	29	BM6UR-290	BM6UL-290	BM6LL-290	BM6LR-290
29+	29.5	29.5	BM6UR-295	BM6UL-295	BM6LL-295	BM6LR-295
30	30	30	BM6UR-300	BM6UL-300	BM6LL-300	BM6LR-300
30+	30.5	30.5	BM6UR-305	BM6UL-305	BM6LL-305	BM6LR-305
31	31	31	BM6UR-310	BM6UL-310	BM6LL-310	BM6LR-310
31+	31.5	31.5	BM6UR-315	BM6UL-315	BM6LL-315	BM6LR-315
32	32	32	BM6UR-320	BM6UL-320	BM6LL-320	BM6LR-320
32+	32.5	32.5	BM6UR-325	BM6UL-325	BM6LL-325	BM6LR-325
33	33	33	BM6UR-330	BM6UL-330	BM6LL-330	BM6LR-330
33+	33.5	33.5	BM6UR-335	BM6UL-335	BM6LL-335	BM6LR-335
34	34	34	BM6UR-340	BM6UL-340	BM6LL-340	BM6LR-340
34+	34.5	34.5	BM6UR-345	BM6UL-345	BM6LL-345	BM6LR-345
35	35	35	BM6UR-350	BM6UL-350	BM6LL-350	BM6LR-350
35+	35.5	35.5	BM6UR-355	BM6UL-355	BM6LL-355	BM6LR-355
36	36	36	BM6UR-360	BM6UL-360	BM6LL-360	BM6LR-360
36+	36.5	36.5	BM6UR-365	BM6UL-365	BM6LL-365	BM6LR-365
37	37	37	BM6UR-370	BM6UL-370	BM6LL-370	BM6LR-370
37+	37.5	37.5	BM6UR-375	BM6UL-375	BM6LL-375	BM6LR-375
38	38	38	BM6UR-380	BM6UL-380	BM6LL-380	BM6LR-380
38+	38.5	38.5	BM6UR-385	BM6UL-385	BM6LL-385	BM6LR-385
39	39	39	BM6UR-390	BM6UL-390	BM6LL-390	BM6LR-390
39+	39.5	39.5	BM6UR-395	BM6UL-395	BM6LL-395	BM6LR-395
40	40	40	BM6UR-400	BM6UL-400	BM6LL-400	BM6LR-400
40+	40.5	40.5	BM6UR-405	BM6UL-405	BM6LL-405	BM6LR-405
41	41	41	BM6UR-410	BM6UL-410	BM6LL-410	BM6LR-410
41+	41.5	41.5	BM6UR-415	BM6UL-415	BM6LL-415	BM6LR-415
42	42	42	BM6UR-420	BM6UL-420	BM6LL-420	BM6LR-420
42+	42.5	42.5	BM6UR-425	BM6UL-425	BM6LL-425	BM6LR-425
43	43	43	BM6UR-430	BM6UL-430	BM6LL-430	BM6LR-430
43+	43.5	43.5	BM6UR-435	BM6UL-435	BM6LL-435	BM6LR-435
44	44	44	BM6UR-440	BM6UL-440	BM6LL-440	BM6LR-440
44+	44.5	44.5	BM6UR-445	BM6UL-445	BM6LL-445	BM6LR-445



Molar Band Prewelded with Tubes

All molar bands can be pre-welded with lingual attachments, and "X" stands for lingual attachments options:

1. Lingual Cleats

A. lingual cleats big



B. lingual cleats small

2. Lingual Sheath

A. Lingual sheath no hook no slot

B. Lingual sheath no hook with slot

C. Lingual sheath mesial with hook no slot

D. Lingual sheath distal with hook no slot

E. Lingual sheath mesial with hook with slot

F. Lingual sheath distal with hook with slot



3. Lingual Button

A. Lingual button round



B. Lingual button rectangular

1st Molar Band Pre-Welded With Single Tube (Individual Pack)

Prescription	Teeth	Torque	Offset	Width (mm)	Item Number (Regular: No Funnel Entrance / Round Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6°	4	BM131-1126-X / BM231-1126-X	BM132-1126-X / BM232-1126-X	BM131-2126-X / BM231-2126-X	BM132-2126-X / BM232-2126-X
	UL6	0°	6°	4	BM131-1116-X / BM231-1116-X	BM132-1116-X / BM232-1116-X	BM131-2116-X / BM231-2116-X	BM132-2116-X / BM232-2116-X
	LL6	0°	6°	4	BM131-1136-X / BM231-1136-X	BM132-1136-X / BM232-1136-X	BM131-2136-X / BM231-2136-X	BM132-2136-X / BM232-2136-X
	LR6	0°	6°	4	BM131-1146-X / BM231-1146-X	BM132-1146-X / BM232-1146-X	BM131-2146-X / BM231-2146-X	BM132-2146-X / BM232-2146-X
Roth	UR6	-14°	10°	4	BM111-1126-X / BM211-1126-X	BM112-1126-X / BM212-1126-X	BM111-2126-X / BM211-2126-X	BM112-2126-X / BM212-2126-X
	UL6	-14°	10°	4	BM111-1116-X / BM211-1116-X	BM112-1116-X / BM212-1116-X	BM111-2116-X / BM211-2116-X	BM112-2116-X / BM212-2116-X
	LL6	-25°	4°	4	BM111-1136-X / BM211-1136-X	BM112-1136-X / BM212-1136-X	BM111-2136-X / BM211-2136-X	BM112-2136-X / BM212-2136-X
	LR6	-25°	4°	4	BM111-1146-X / BM211-1146-X	BM112-1146-X / BM212-1146-X	BM111-2146-X / BM211-2146-X	BM112-2146-X / BM212-2146-X
MBT	UR6	-14°	10°	4	BM121-1126-X / BM221-1126-X	BM122-1126-X / BM222-1126-X	BM121-2126-X / BM221-2126-X	BM122-2126-X / BM222-2126-X
	UL6	-14°	10°	4	BM121-1116-X / BM221-1116-X	BM122-1116-X / BM222-1116-X	BM121-2116-X / BM221-2116-X	BM122-2116-X / BM222-2116-X
	LL6	-20°	0°	4	BM121-1136-X / BM221-1136-X	BM122-1136-X / BM222-1136-X	BM121-2136-X / BM221-2136-X	BM122-2136-X / BM222-2136-X
	LR6	-20°	0°	4	BM121-1146-X / BM221-1146-X	BM122-1146-X / BM222-1146-X	BM121-2146-X / BM221-2146-X	BM122-2146-X / BM222-2146-X



1st Molar Band Pre-Welded With Single Tube (One Set 4 pcs / Pack)

Prescription	Teeth	Torque	Offset	Width (mm)	Item Number (Regular: No Funnel Entrance / Round Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6°	4	BM131-11B-X / BM231-11B-X	BM132-11B-X / BM232-11B-X	BM131-21B-X / BM231-21B-X	BM132-21B-X / BM232-21B-X
	UL6	0°	6°	4				
	LL6	0°	6°	4				
	LR6	0°	6°	4				
Roth	UR6	-14°	10°	4	BM111-11B-X / BM211-11B-X	BM112-11B-X / BM212-11B-X	BM111-21B-X / BM211-21B-X	BM112-21B-X / BM212-21B-X
	UL6	-14°	10°	4				
	LL6	-25°	4°	4				
	LR6	-25°	4°	4				
MBT	UR6	-14°	10°	4	BM121-11B-X / BM221-11B-X	BM122-11B-X / BM222-11B-X	BM121-21B-X / BM221-21B-X	BM122-21B-X / BM222-21B-X
	UL6	-14°	10°	4				
	LL6	-20°	0°	4				
	LR6	-20°	0°	4				



1st Molar Band Pre-Welded With Double Combination Tube (Individual Pack)

Prescription	Teeth	Torque	Offset	Width (mm)	Item Number (Regular: No Funnel Entrance / Round Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6°	4	BM131-1226-X / BM231-1226-X	BM132-1226-X / BM232-1226-X	BM131-2226-X / BM231-2226-X	BM132-2226-X / BM232-2226-X
	UL6	0°	6°	4	BM131-1216-X / BM231-1216-X	BM132-1216-X / BM232-1216-X	BM131-2216-X / BM231-2216-X	BM132-2216-X / BM232-2216-X
	LL6	0°	6°	4	BM131-1236-X / BM231-1236-X	BM132-1236-X / BM232-1236-X	BM131-2236-X / BM231-2236-X	BM132-2236-X / BM232-2236-X
	LR6	0°	6°	4	BM131-1246-X / BM231-1246-X	BM132-1246-X / BM232-1246-X	BM131-2246-X / BM231-2246-X	BM132-2246-X / BM232-2246-X
Roth	UR6	-14°	10°	4	BM111-1226-X / BM211-1226-X	BM112-1226-X / BM212-1226-X	BM111-2226-X / BM211-2226-X	BM112-2226-X / BM212-2226-X
	UL6	-14°	10°	4	BM111-1216-X / BM211-1216-X	BM112-1216-X / BM212-1216-X	BM111-2216-X / BM211-2216-X	BM112-2216-X / BM212-2216-X
MBT	LL6	-20°	0°	4	BM121-1226-X / BM221-1226-X	BM122-1226-X / BM222-1226-X	BM121-2226-X / BM221-2226-X	BM122-2226-X / BM222-2226-X
	LR6	-20°	0°	4	BM121-1216-X / BM221-1216-X	BM122-1216-X / BM222-1216-X	BM121-2216-X / BM221-2216-X	BM122-2216-X / BM222-2216-X



1st Molar Band Pre-Welded With Double Combination Tube (One Set 4 pcs / Pack)

Prescription	Teeth	Torque	Offset	Width (mm)	Item Number (Regular: No Funnel Entrance / Round Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6°	4	BM131-12B-X / BM231-12B-X	BM132-12B-X / BM232-12B-X	BM131-22B-X / BM231-22B-X	BM132-22B-X / BM232-22B-X
	UL6	0°	6°	4				
	LL6	0°	6°	4				
	LR6	0°	6°	4				
Roth	UR6	-14°	10°	4	BM111-12B-X / BM211-12B-X	BM112-12B-X / BM212-12B-X	BM111-22B-X / BM211-22B-X	BM112-22B-X / BM212-22B-X
	UL6	-14°	10°	4				
	LL6	-25°	4°	4				
	LR6	-25°	4°	4				
MBT	UR6	-14°	10°	4	BM121-12B-X / BM221-12B-X	BM122-12B-X / BM222-12B-X	BM121-22B-X / BM221-22B-X	BM122-22B-X / BM222-22B-X
	UL6	-14°	10°	4				
	LL6	-20°	0°	4				
	LR6	-20°	0°	4				



1st Molar Band Pre-Welded With Double Rectangular Tube (Individual Pack)

Prescription	Teeth	Torque	Offset	Width (mm)	Item Number (Regular: No Funnel Entrance / Round Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6°	4	BM131-1226-X / BM231-1226-X	BM132-1226-X / BM232-1226-X	BM131-2226-X / BM231-2226-X	BM132-2226-X / BM232-2226-X
	UL6	0°	6°	4	BM131-1216-X / BM231-1216-X	BM132-1216-X / BM232-1216-X	BM131-2216-X / BM231-2216-X	BM132-2216-X / BM232-2216-X
	LL6	0°	6°	4	BM131-1236-X / BM231-1236-X	BM132-1236-X / BM232-1236-X	BM131-2236-X / BM231-2236-X	BM132-2236-X / BM232-2236-X
	LR6	0°	6°	4	BM131-1246-X / BM231-1246-X	BM132-1246-X / BM232-1246-X	BM131-2246-X / BM231-2246-X	BM132-2246-X / BM232-2246-X
Roth	UR6	-14°	10°	4	BM111-1226-X / BM211-1226-X	BM112-1226-X / BM212-1226-X	BM111-2226-X / BM211-2226-X	BM112-2226-X / BM212-2226-X
	UL6	-14°	10°	4	BM111-1216-X / BM211-1216-X	BM112-1216-X / BM212-1216-X	BM111-2216-X / BM211-2216-X	BM112-2216-X / BM212-2216-X
	LL6	-25°	4°	4	BM111-1236-X / BM211-1236-X	BM112-1236-X / BM212-1236-X	BM111-2236-X / BM211-2236-X	BM112-2236-X / BM212-2236-X
	LR6	-25°	4°	4	BM111-1246-X / BM211-1246-X	BM112-1246-X / BM212-1246-X	BM111-2246-X / BM211-2246-X	BM112-2246-X / BM212-2246-X
MBT	UR6	-14°	10°	4	BM121-1226-X / BM221-1226-X	BM122-1226-X / BM222-1226-X	BM121-2226-X / BM221-2226-X	BM122-2226-X / BM222-2226-X
	UL6	-14°	10°	4	BM121-1216-X / BM221-1216-X	BM122-1216-X / BM222-1216-X	BM121-2216-X / BM221-2216-X	BM122-2216-X / BM222-2216-X
	LL6	-20°	0°	4	BM121-1236-X / BM221-1236-X	BM122-1236-X / BM222-1236-X	BM121-2236-X / BM221-2236-X	BM122-2236-X / BM222-2236-X
	LR6	-20°	0°	4	BM121-1246-X / BM221-1246-X	BM122-1246-X / BM222-1246-X	BM121-2246-X / BM221-2246-X	BM122-2246-X / BM222-2246-X



1st Molar Band Pre-Welded With Double Rectangular Tube (One Set 4 pcs / Pack)

Prescription	Teeth	Torque	Offset	Width (mm)	Item Number (Regular: No Funnel Entrance / Round Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6°	4	BM131-12B-X / BM231-12B-X	BM132-12B-X / BM232-12B-X	BM131-22B-X / BM231-22B-X	BM132-22B-X / BM232-22B-X
	UL6	0°	6°	4				
	LL6	0°	6°	4				
	LR6	0°	6°	4				
Roth	UR6	-14°	10°	4	BM111-12B-X / BM211-12B-X	BM112-12B-X / BM212-12B-X	BM111-22B-X / BM211-22B-X	BM112-22B-X / BM212-22B-X
	UL6	-14°	10°	4				
	LL6	-25°	4°	4				
	LR6	-25°	4°	4				
MBT	UR6	-14°	10°	4	BM121-12B-X / BM221-12B-X	BM122-12B-X / BM222-12B-X	BM121-22B-X / BM221-22B-X	BM122-22B-X / BM222-22B-X
	UL6	-14°	10°	4				
	LL6	-20°	0°	4				
	LR6	-20°	0°	4				



1st Molar Band Pre-Welded With Triple Tube (Individual Pack)

Prescription	Teeth	Torque	Offset	Width (mm)	Item Number (Regular: No Funnel Entrance / Round Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6°	4	BM131-1326-X / BM231-1326-X	BM132-1326-X / BM232-1326-X	BM131-2326-X / BM231-2326-X	BM132-2326-X / BM232-2326-X
	UL6	0°	6°	4	BM131-1316-X / BM231-1316-X	BM132-1316-X / BM232-1316-X	BM131-2316-X / BM231-2316-X	BM132-2316-X / BM232-2316-X
Roth	UR6	-14°	10°	4	BM111-1326-X / BM211-1326-X	BM112-1326-X / BM212-1326-X	BM111-2326-X / BM211-2326-X	BM112-2326-X / BM212-2326-X
	UL6	-14°	10°	4	BM111-1316-X / BM211-1316-X	BM112-1316-X / BM212-1316-X	BM111-2316-X / BM211-2316-X	BM112-2316-X / BM212-2316-X
MBT	LL6	-20°	0°	4	BM121-1326-X / BM221-1326-X	BM122-1326-X / BM222-1326-X	BM121-2326-X / BM221-2326-X	BM122-2326-X / BM222-2326-X
	LR6	-20°	0°	4	BM121-1316-X / BM221-1316-X	BM122-1316-X / BM222-1316-X	BM121-2316-X / BM221-2316-X	BM122-2316-X / BM222-2316-X



1st Molar Band Pre-Welded With Triple Tube (One Set 4 pcs / Pack)

Prescription	Teeth	Torque	Offset	Width (mm)	Item Number (Regular: No Funnel Entrance / Round Funnel Entrance)			
					Non-Convertible		Convertible	
					0.018	0.022	0.018	0.022
Egewise	UR6	0°	6°	4	BM131-13B-X / BM231-13B-X	BM132-13B-X / BM232-13B-X	BM131-23B-X / BM231-23B-X	BM132-23B-X / BM232-23B-X
	UL6	0°	6°	4				
Roth	UR6	-14°	10°	4	BM111-13B-X / BM211-13B-X	BM112-13B-X / BM212-13B-X	BM111-23B-X / BM211-23B-X	BM112-23B-X / BM212-23B-X
	UL6	-14°	10°	4				
MBT	LL6	-20°	0°	4	BM121-13B-X / BM221-13B-X	BM122-13B-X / BM222-13B-X	BM121-23B-X / BM221-23B-X	BM122-23B-X / BM222-23B-X
	LR6	-20°	0°	4				



Orthodontic Accessories



Product Description	Size / Style	Item Number	Packing
Crimpable Hooks	Round	CHS-11	10 pcs / pcak
	Rectangular	CHS-12	10 pcs / pcak
Crimpable Hooks Long	Round	CHSL-11	10 pcs / pcak
	Rectangular	CHSL-12	10 pcs / pcak
Crimpable Stop	Round	CSS-11	10 pcs / pcak
	Rectangular	CSS-12	10 pcs / pcak
Long Curved Crimpable Hook	Left	LCCHS-L	10 pcs / pcak
	Right	LCCHS-R	10 pcs / pcak
Lingual Buttons	Bondable	LBS-B	10 pcs / pcak
	Weldable	LBS-W	10 pcs / pcak
Lingual Cleats	Bondable	LCT-B	10 pcs / pcak
	Weldable	LCT-W	10 pcs / pcak
Lingual Sheath	no hook no slot	LSH-A	10 pcs / pcak
	no hook with slot	LSH-B	10 pcs / pcak
	mesial with hook no slot	LSH-C	10 pcs / pcak
	distal with hook no slot	LSH-D	10 pcs / pcak
	mesial with hook with slot	LSH-E	10 pcs / pcak
	distal with hook with slot	LSH-F	10 pcs / pcak
Lingual Retainers	38 - Upper	LRS-U38	10 pcs / pcak
	40 - Upper	LRS-U40	
	42 - Upper	LRS-U42	
	45 - Upper	LRS-U45	
	47 - Upper	LRS-U47	10 pcs / pcak
	29 - Lower	LRS-L29	
	31 - Lower	LRS-L31	
	33 - Lower	LRS-L33	
35 - Lower	LRS-L35		
37 - Lower	LRS-L37		
Standard Face Bow	/	FB01	10 pcs / pcak
Loop Style Face Bow	83mm	FB02-1	10 pcs / pcak
	90mm	FB02-2	10 pcs / pcak
	97mm	FB02-3	10 pcs / pcak
	104mm	FB02-4	10 pcs / pcak
	111mm	FB02-5	10 pcs / pcak

